

Genetics 284 Laboratory Schedule (Tentative)

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	1/24	Intro to Flies		<p>Monday night clear vials of wildtype Canton-S, put at 18C for virgins in the morning.</p> <p>12 vials of wildtype Canton-S 12 vials of mutant flies, mapping stocks</p>
T	1/31	Intro to Worms		<p>12 plates of N2 normal worms 1 plate of each mutation</p> <p>Alcohol lamps (12) Picks (12 at least - ask Dave)</p>

Day	Date	Lab Day	Outside class	Lab Prep Needs
M	2/6		<p>Start overnight of RNAi plasmid containing bacteria</p> <p>pL4440 empty pL4440 + lsy-2 pL4440 + choice</p> <p>37C incubator</p>	<p>Plates of bacteria with:</p> <p>pL4440 empty pL4440 + lsy-2 (chemo assay)</p> <p>Groups will choose one of the following for phenotypic analysis:</p> <p>pL4440 + bli-1 pL4440 + dpy-5 pL4440 + unc-22 (from Chase Lab)</p> <p>10ml per overnight of sterile LB + 12.5ug/ml tetracycline</p> <p>Adding 10ul of 50mg/ml ampicillin stock</p> <p>15ml conical tubes for overnight, 3 per group, so 36 tubes = 360ml of broth minimum</p> <p>Sterile toothpicks</p>

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	2/7	<p>Set up PCR to verify plasmid inserts</p> <p>Isolate plasmid DNA from bacteria</p>	Begin collecting virgin females from the P-element lines this week	<p>PCR reagents:</p> <p>Invitrogen T7 Promoter Primer Catalog Number N560-02 (forward primer, 20uM)</p> <p>50 ul of ultrapure sterile water</p> <p>PCR buffer 50mM MgCl₂ 10mM dNTP mix Taq Polymerase (5U/microliter)</p> <p>PCR engines Parameters: 94C 5min 94C 30sec 55C 30sec 72C 3min Repeat 30 72C 6min 4C HOLD</p> <p>Solutions for DNA isolation:</p> <p>Solution 1 (~ 7.2 mls min) 50 mM glucose 10 mM EDTA 25 mM Tris, pH 8</p> <p>Solution 2: (~14.4 mls min) 0.2 M NaOH 1% SDS (sodium dodecyl sulfate)</p> <p>Solution 3: (~10.8mls min) 2.7M Potassium Acetate 6.6M Acetic Acid</p> <p>-----</p> <p>100% ethanol 80% ethanol sterile water Speed Vac??</p>

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	2/14	<p>Gel analysis of PCR products, verification of inserts</p> <p>Transformation of HTT115DE3 feeding lines, plates to 37C</p> <p>-----</p> <p>Set up Cross 1 p-element virgins x males for mapping stock (18C incubator)</p>		<p>2% agarose 1 x TAE mini gels to check PCR products (EtBR)</p> <p>loading dye for samples?</p> <p>Need competent HTT115DE3 E.coli made by this date (protocol by Inoue, see at end of this document)</p> <p>LB broth LB + amp + tet plates</p> <p>-----</p> <p>Dissecting Scopes out, fly morgue, anesthesitizers, Fly Nap, brushes, clean vials with food</p> <p>Multiple marker chromosome mapping stocks to get males from.</p> <p>18C incubator</p>
W	2/15		<p>Check transformation and control plates</p> <p>Save transformation plates, notify TA/instructor if control plate is positive also!</p>	
M	2/20		<p>Overnight cultures of transformed bacteria</p>	<p>LB + tet broth amp stock sterile toothpicks</p>

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	2/21	<p>Induce HTT115DE3 bacteria to produce RNAi by IPTG</p> <p>Seed plates with induced bacteria for feeding plates, allow plates to dry</p> <p>-----</p> <p>Clear parent flies from Cross 1 vials put back in 18C incubator</p>		<p>Overnight cultures will need to be subcultured by Prep Staff in the morning in order to be in log phase by lab time for induction</p> <p>Diluting cultures 1:10 (500ul culture into 4.5ml of LB+tet+amp) put in 37C to grow 3-4 hours until lab time.</p> <p>0.5M IPTG stock</p> <p>NGM plates with 0.4mM IPTG, 50ug/ml ampicillin, and 12.5ug/mL tetracycline.</p>
Fri	2/24	Put N2 and rrf-3 worms on the feeding plates to feed them RNAi		<p>L4s of N2 worms rrf-3 worms</p> <p>23C incubator</p>

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	2/28	<p>Collect phenotypic data on gene of interest</p> <p>Wash worms for chemotaxis assay from feeding plates, resuspend and put on chemotaxis plates</p> <p>Chemotaxis assay</p>		<p>Day before chemotaxis assay plates set up. (Prep staff or students??)</p> <p>Experimental: NGM plates no food, 10ul of 2.5M NaCl on one side of plate, 10ul of double deionized water for control, mark bottom of plate with location of each</p> <p>Control: have water on both sides marked</p> <p>ice cold sterile water</p> <p>0.25M sodium azide 2.5M NaCl</p>

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	3/6	<p>Isolate Fly DNA from P-element line</p> <p>Digest DNA to cut out the P-element</p> <p>Ligate o/n at 4C</p>		<p>Isolating DNA from flies: Need bottles of P-element lines by here, 30 flies per group</p> <p>Solution A: (~6 ml) 100mM Tris-HCL pH 7.5 100mM EDTA</p> <p>LiCl/KAc solution: (~12ml) 1 part 5M KAc stock 2.5 parts 6M LiCl stock</p> <p>65C heat block/water bath 37C heat block/water bath</p> <p>isopropanol 70% ethanol TE for resuspension</p> <p>Digest: 10X RE buffer 100um/ml RNase ddH2O Sau3A I, HinP1 I or Msp I (Kenny? which one for P[LacW])</p> <p>Ligation: 10X Ligation buffer (+ATP) ddH2O Ligase (2 Weiss units)</p>
W	3/7		Precipitate DNA and resuspend in 150ul TE	<p>Ethanol ice Speed Vac TE for resuspension</p> <p>(Students do?? TA does??)</p>

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	3/13	<p>Run PCR to amplify flanking sequence from P-element primers</p> <p>Mini gel check for amplification, put in 4C</p> <p>Send out sequencing here?</p>		<p>PCR on P-element DNA: 2mM each dNTP 10uM forward primer 10uM reverse primer 10X Pharmacia Taq buffer ddH2O 2 U Taq</p> <p>PCR parameters: 95C 5min</p> <p>95C 30sec 60 or 55C 1 min 68C 2min</p> <p>Repeat 35X 72C 10 min 4C HOLD</p> <p>Gel Check PCR products: 1.5% agarose mini gel with EtBr markers??</p> <p>Who sets up sequencing reaction to go to GeneWiz?? Students or TA??</p>
T	3/20	SPRING BREAK	SPRING BREAK	SPRING BREAK
T	3/27	<p>Bioinformatics Session, use sequence to locate insertion site</p> <p>Sequence back here?</p>		

Day	Date	Lab Day	Outside class	Lab Prep Needs
T	4/3	Score flies from the mapping cross, do three point cross and determine recombinant map position		fly sorting dishes! morgue
T	4/10	LaZ staining of P-element lines??		
T	4/17	HOLIDAY	HOLIDAY	HOLIDAY
T	4/24	Presentations		
T	5/1	Presentations		